



## Benedict's Qualitative Reagent

R002

### Intended use

Benedict's reagent is used to test for the presence of reducing sugars.

### Composition\*\*

#### Ingredients

Copper sulphate	17.300 gm
Sodium carbonate	100.000 gm
Sodium citrate	173.000 gm
Distilled water	1000.000 ml

\*\*Formula adjusted, standardized to suit performance parameters

### Directions

For detection of sugar in Urine:

Add 5 ml of Benedict's qualitative reagent in a test tube. Add 8 drops (0.5 ml) of urine. Boil over a flame (or in a boiling water bath) for 5-10 minutes. Cool under tap water. The contents of the tube becomes turbid due to a precipitate, which may range from green to brick red in colour, depending on the amount of sugar present in the urine. If no sugar is present, the solution will remain clear or show a faint turbidity.

### Principle And Interpretation

Benedict's reagent is used to test for the presence of glucose in urine. Once a reducing sugar is detected in urine, further tests have to be undergone in order to ascertain which sugar is present. The copper sulphate in Benedict's solution reacts with reducing sugars and the cupric ions to cuprous ions, these are precipitated as red copper oxide, which is insoluble in water (1). Alkaline medium is provided to the reaction by sodium carbonate present in the reagent. The original colour of Benedict's reagent is blue. It changes to green, yellow, orange or red, according to the concentration of glucose present in urine.

### Type of specimen

Clinical samples

### Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (2,3).

After use, contaminated materials must be sterilized by autoclaving before discarding.

### Warning and Precautions :

In Vitro diagnostic Use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets

### Performance and Evaluation

Performance of the product is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

### Quality Control

#### Appearance

Blue coloured solution.

#### Clarity

Clear to very slightly opalescent solution.

#### Test

Procedure: Add 5 ml of Benedict's qualitative reagent in a test tube. Add 8 drops (0.5 ml) of urine. Boil over a flame for 3 min. Cool. The contents of tube become turbid due to precipitate which may range from green to brick red in colour depending on the amount of sugar present in the urine. If no sugar is present, the solution will remain clear or show a faint turbidity.

**Result**

Colour of Mixture	Approximate amount of glucose	Conclusion (Sugar)
Blue	Nil	Absent
Green	0.5%	Present, trace
Greenish brown	1.0%	Present, + to ++
Yellow	1.5%	Present, +++
Brick red	2.0% or more	Present, ++++

**Storage and Shelf Life**

Store between 10-30°C in tightly closed container and away from bright light. Use before expiry date on label. On opening, product should be properly stored in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use.

**Disposal**

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (2,3).

**Reference**

1. Benedict, S.R. "A Reagent for the detection of Reducing Sugars", J. Biol. Chem. 5(6):485-487.
2. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2<sup>nd</sup> Edition.
3. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S. and Warnock, D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.


Revision : 02 / 2019



In vitro diagnostic  
medical device



CE  
Marking



Storage  
temperature



Do not use if  
package is  
damaged



HiMedia Laboratories Pvt.  
Limited, B /4-6 , MIDC,  
Dindori, Nashik MH  
[www.himedialabs.com](http://www.himedialabs.com)



CE Partner 4U ,Esdoornlaan  
13, 3951 DB Maarn The  
Netherlands,  
[www.cepartner4u.eu](http://www.cepartner4u.eu)

**Disclaimer :**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.