



## DMACA Reagent

R035

### Intended Use

DMACA Reagent is used as a histological dye to detect indole from Tryptophan.

### Composition\*\*

#### Ingredients

p-Dimethylaminocinnamaldehyde	1.000 gm
Hydrochloric acid (concentrated)	1.000 ml
Distilled water	99.000 ml

\*\*Formula adjusted, standardized to suit performance parameters

### Directions

With a platinum wire loop, plastic loop or glass rod spread suspected colony from HiCrome UTI HiVeg Agar (MV1353) or HiCrome UTI HiVeg Agar, Modified (MV1418) plate on filter paper dipped in DMACA Reagent (R035). Observe for appearance of bluish-purple colour within 10-30 seconds.

### Principle And Interpretation

DMACA reagent is used to determine the ability of an organism to split indole from the tryptophan molecule. In the presence of oxygen, some bacteria are able to split tryptophan into indole and alpha-aminopropionic acid. The presence of indole can be detected by the addition of DMACA (p-Dimethylaminocinnamaldehyde) reagent indicated by formation of bluish-purple colour.

### Type of specimen

The specimen is any isolated colony on primary or subculture plates

### Specimen Collection and Handling

or clinical samples follow appropriate techniques for handling specimens as per established guidelines

or food and dairy samples follow appropriate techniques for sample collection and processing as per guidelines 3.5

3 or watersamples follow appropriate techniques for sample collection processing as per guidelines and local standards

After use contaminated materials must be sterilized by autoclaving before discarding

### Warning and Precautions

In vitro diagnostic use only Read the label before opening the container Wear protective gloves protecting clothing eye protection face protection Follow good microbiological lab practices while handling specimens and cultures Standard

### Limitations

1. The tube test is a more sensitive method of detecting indole than the spot test.
2. When performing a spot test, Kovacs Indole Reagent may be used as a substitute for the spot test reagent. However, Kovacs Indole Reagent, when used as the spot test reagent, is less sensitive in detecting indole than the Indole Spot Reagent (DMACA).
3. Media containing glucose should not be used for indole testing due to the formation of acid end products which have been shown to reduce indole production.
4. Indole-positive colonies have been reported to cause adjacent indole-negative colonies to appear false-positive due to diffusion of indole into the media. To avoid false-positives, select colonies of different morphologies that are separated by at least 5mm for indole testing.

---

## Quality Control

### Appearance

Yellow to brown coloured solution.

### Clarity

Clear to slightly hazy solution.

### Cultural Response

Biochemical identification is carried out by transferring suspected colony from HiCrome UTI Agar (M1353) or HiCrome UTI Agar, Modified (M1418) plate on filter paper and adding 1-2 drops of DMACA Reagent (R035).

Organism	Growth	DMACA
<i>Escherichia coli</i> ATCC 25922	Luxuriant	Positive (blue-purple colouration around colony)
<i>Pseudomonas aeruginosa</i> ATCC 27853	Luxuriant	Negative(No colour change around colony)

## Storage and Shelf Life

Store below 30°C in tightly closed container and away from bright light. Use before expiry date on label.

## Reference

1. MacFaddin J. F., 2000, Biochemical Tests for Identification of Medical Bacteria, 3rd Ed., Williams and Wilkins, Baltimore.
2. Bogs, J.; Jaffe, F. W.; Takos, A. M.; Walker, A. R.; Robinson, S. P. (2007). "The Grapevine Transcription Factor VvMYBPA1 Regulates Proanthocyanidin Synthesis during Fruit Development". *Plant Physiology* 143 (3): 1347–1361.
3. Li, Y. G.; Tanner, G.; Larkin, P. (1996). "The DMACA-HCl Protocol and the Threshold Proanthocyanidin Content for Bloat Safety in Forage Legumes". *Journal of the Science of Food and Agriculture* 70: 89–101.
4. A new colourimetric assay for flavonoids in pilsner beers. Jan A. Delcour and Didier Janssens de Varebeke, *Journal of the Institute of Brewing*, January–February 1985, Volume 91, Issue 1, pages 37–40.
5. Meudt, W. J.; Gaines, T. P. (1967). "Studies on the Oxidation of Indole-3-Acetic Acid by Peroxidase Enzymes. I. Colorimetric Determination of Indole-3-Acetic Acid Oxidation Products". *Plant Physiology* 42 (10): 1395–9.

Revision : 02 / 2016



### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.