

Trypsin 1:250 powder porcine

Gamma Irradiated

Source: Porcine pancreas

Activity: 1000 - 1500 BAEE units/mg

Cell Culture Tested

Product Code: TC245G

Molecular weight: 23.4 kDA

CAS No: 9002-07-7

EC No: 3.4.21.4

Extinction Co-efficient: $E_{1\%}^{1\text{cm}} = 14.3(280\text{nm})^2$

Trypsin is a serine protease derived from pancreas. It is a single chain polypeptide of 223 amino acid residues with substrate specificity based on positively charged Lysine and Arginine side chains. Trypsin predominantly cleaves peptide chains at the carboxyl side of Lysine and Arginine, except when either is followed by Proline.

Trypsin is produced from Trypsinogen by removal of a terminal hexapeptide to yield a single chain native form of trypsin called β -Trypsin. Subsequent autolysis of β -Trypsin results in α -Trypsin having two peptide chains bound by disulphide bonds.

TC245G is Gamma irradiated trypsin. Gamma irradiation process provides greater assurance that any existing low level of microorganisms will be inactivated or reduced and the risk associated with animal-derived components are minimized.

Activity

The optimum pH for trypsin activity is 7.0-9.0 and optimum temperature is 37°C.

One BAEE unit will produce a r A253nm of 0.001 per minute with BAEE as substrate at pH 7.6 at 25°C in a reaction volume of 3.2ml (1cm light path).

One TAME unit hydrolyzes 1umole of p-toluene-sulfonyl-L-arginine methyl ester (TAME) per minute at 25°C, pH 8.2 in the presence of 0.001M calcium ion.

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ Publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

One USP trypsin unit is the activity causing a change in absorbance of 0.003 per minute under the conditions specified.

Activity Conversion: 1 TAME unit = 19.2 USP or NF units= 57.5 BAEE Units.

Quality Control:

Appearance

White to offwhite powder.

Solubility

Clear colorless solution at 2.5gm in 100ml of Hank's balanced salt solution .

Loss on drying

NMT 5%

Porcine parvovirus

Not detected

Chymotrypsin assay

70.00 -140.00

Cell Culture Test

Passes

Storage and Shelf Life:

Trypsin powder is stable at 2-8°C for 36 months.

Trypsin after reconstitution should be stored at -20°C.

Autolysis of Trypsin may also be prevented by presence of Ca^{+2} ions.