



Grams Stain-Kit

K001

Intended Use

Grams Stain Kit is used for differentiation of bacteria on the basis of their gram nature.

Composition**

Ingredients

Gram's Crystal Violet (S012)(Solution A)

Crystal Violet	2.000 gm
Ethyl alcohol,95%	20.000 ml

Gram's Crystal Violet (S012)(Solution B)

Ammonium oxalate	0.800 gm
Distilled Water	80.000 ml

Solution A and B mixed .Stored for 24 hours before use.The resulting stain is stable.

Gram's Decolourizer(S032)

Ethyl alcohol, 95%	50.0 ml
Acetone	50.0 ml

Gram's Iodine(S013)

Iodine	1.000 gm
Potassium iodide	2.000 gm
Distilled water	300.000 ml

Safranin,0.5% w/v(S027)

Safranin O	0.500 gm
Ethyl alcohol, 95%	100.000 ml

**Formula adjusted, standardized to suit performance parameters

Directions

- 1)Prepare a thin smear on clear, dry glass slide.
- 2)Allow it to air dry and fix by gentle heat.
- 3)Flood with Gram's Crystal Violet (S012) for 1 minute. (If over staining results in improper decolourization of known gram-negative organisms,use less crystal violet).
- 4)Drain the stain.
- 5)Flood the smear with Gram's Iodine (S013). Allow it to remain for 1 minute.
- 6)Decolourize with Gram's Decolourizer (S032) until the blue dye no longer flows from the smear.
- 7)Wash with tap water.
- 8)Counter stain with 0.5% w/v Safranin (S027). Allow it to remain for 1 minute.
- 9)Wash with water.
- 10)Allow the slide to air dry or blot dry between sheets of clean bibulous paper and examine under oil immersion objective.

Principle And Interpretation

The Gram stain is a differential staining technique most widely applied in all microbiology disciplines laboratories. It is one of the most important criteria in any identification scheme for all types of bacterial isolates. Different mechanisms have been proposed to explain the gram reaction. There are many physiological differences between gram-positive and gram-negative cell walls (1). Ever since Christian Gram has discovered Gram staining, this process has been extensively investigated and redefined. In practice,a thin smear of bacterial cells is stained with crystal violet, then treated with an iodine containing mordant to increase the binding of primary stain (2). A decolourizing solution of alcohol or acetone is used to remove the crystal violet from cells which bind it weakly and then the counterstain (like safranin) is used to provide a colour contrast in those cells that are decolourized.

Gram-positive bacteria have a thick mesh-like cell wall made of peptidoglycan (50–90% of cell envelope), and as a result are stained purple by crystal violet, whereas gram-negative bacteria have a thinner layer (10% of cell envelope), so do not retain the purple stain and are counter-stained pink by safranin. In a properly stained smear by gram staining procedure, the gram-positive bacteria appear blue to purple and gram negative cells appear pink to red.

Type of specimen

Clinical samples - Blood, urine, CSF, pus, wounds, lesions, body tissues, sputum etc. ; food & dairy samples ; Water samples

Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (4,5).

For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards(2).

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions :

In Vitro diagnostic Use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets

Limitations :

1. Growth media must contain an adequate amount of tryptophan. Do not use Mueller- Hinton Agar for test, because tryptophan is destroyed during the acid hydrolysis of casein.
2. Do not use media that contain dyes (e.g., EMB, MAC).
3. Do not use medium with a nitrate disc/strip to perform the indole test, as nitrate can interfere with indole test by including false positive results.

Performance and Evaluation

Performance of the product is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Microscopic examination

Gram staining is carried out and observed under oil immersion lens.

Results

Gram-positive organisms : Violet coloured

Gram-negative organisms : Pinkish red coloured

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

Storage and Shelf Life

Store below 30°C in tightly closed container and away from bright light. Use before expiry date on label. On opening, product should be properly stored in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques (6,7).

Reference

- 1.Lamanna and Mallette, 1965, Basic Bacteriology, 3rd ed., Williams and Wilkins Co.,Baltimore.
- 2.Salton, 1964, The Bacterial Cell Wall, Elsevier, Amsterdam.


Revision : 01 / 2017



In vitro diagnostic medical device



CE Marking



Storage temperature



Do not use if package is damaged



HiMedia Laboratories Pvt. Limited,
B /4-6 , MIDC, Dindori, Nashik MH

www.himedialabs.com



CE Partner 4U ,Esdoornlaan 13, 3951
DB Maarn The Netherlands,
www.cepartner4u.eu



Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory,diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.