

## MBRE016

## Xba I

### Components

Reagents provided	MBRE016			
	250 Units	500 Units	1000 Units	5000 Units
Xba I	12.5 µl	25 µl	50 µl	250 µl
10X HiBuffer H5	312.5 µl	625 µl	1.25 ml	6.25 ml
10X HiBuffer DB	312.5 µl	625 µl	1.25 ml	6.25 ml
Diluent E Buffer	162.5 µl	325 µl	650µl	3.25 ml

**NOTE: BSA included in all Reaction Buffer**

**Source:** A *E.coli* strain that carries the Xba I gene from *Xanthomonas badrii*

### Recognition Sequence:



**Concentration:** 20 U/µl

### Unit Definition:

1 u is defined as the amount of enzyme that is required to digest 1µg of DNA in 1 hour at 37°C in 50µl of assay buffer.

Enzyme	Optimum reaction temperature (°C)	Thermal Inactivation (°C)	% activity of Buffers				
			H1	H2	H3	H4	H5
Xba I	37	65	10	75	75	10	100

### Reaction Buffer:

#### 10X HiBuffer H5:

30mM Tris –acetate (pH 7.9 at 30°C), 10 mM Mg-acetate, 60 mM K-acetate and 100 µg /ml BSA.

**NOTE:** Incubate at 37°C

#### Storage Buffer:

10mM Tris-HCl (pH 7.5), 300mM NaCl, 0.1mM EDTA, 7mM 2-mercaptoethanol, 200µg/ml BSA and 50% glycerol. Store at –20°C .

**NOTE: 10X HiBuffer DB is provided for double digestion.**

### Quality Control Assays:

#### Ligation / Recutting Assay :

After 20 -fold overdigestion with Xba I, more than 90% of the DNA fragments can be ligated and recut.

#### Over digestion assay :

An unaltered banding pattern was observed after 1µg of DNA was digested with 40U of Xba I for 16

hours at 37°C.

**Example of Digestion conditions:**

- Enzyme concentration : 1 Unit
- Lambda DNA (0.3 µg /µl) : 3.33 µl (1 µg DNA)
- 10X HiBuffer H5 : 5 µl
- Nuclease free water : upto 50 µl

**Note:**

- Total reaction volume is dependent on the experiment
- The amount of enzyme to be used is dependent on the DNA template
- For plasmid DNA, 5-10X more enzyme is required
- High enzyme concentration may result in **Star activity**

**Storage conditions:** Xba I should be stored at -20°C.

MBRE016\_/0611

MBRE016-02