



o-Toluidine Reagent

R024

It is use for estimation of blood glucose.

Composition**

Ingredients

Thiourea	1.50 gm
Glacial acetic acid	940.0 ml
O-Toluidine ,pure	60.0ml

**Formula adjusted, standardized to suit performance parameters

Directions

Prepare one Test, one Standard and one Blank for each assay series and dispense as follows :

	Test	Standard	Blank
1) O-toluidine reagent, (ml)	5.0	5.0	5.0
2) Glucose (test) of unknown concentration, (mg%)	0.05	-	-
3) Glucose std: 100 (mg %)	-	0.05	-
4) Distilled water, (ml)	0.05	0.05	0.05

Mix thoroughly, and place the tubes in boiling waterbath for exactly 10 minutes. By using tap water, cool the tubes to room temperature. Measure the optical densities of test and standard against blank at 640 nm (red filter 620 - 660 nm).

Principle And Interpretation

A solution of O-toluidine in glacial acetic acid when heated with glucose produces a green coloured complex with an absorption maximum at about 630 nm. The aldehyde group of the glucose apparently condenses with the reagent to form a glucosylamine and a Schiff base, which is probably the coloured product.

Quality Control

Appearance

Light yellow liquid becoming reddish brown on exposure to air and light.

Clarity

Clear without any precipitate

Test

Determination of unknown glucose concentration is carried out using O- toluidine reagent for test sample. The optical densities of test and standard is measured against blank at 640 nm (red filter 620 - 660 nm).

Calculation

Conc. of sample (test),mg% = Absorbance of test/Absorbance of standard * Conc. of standard

Storage and Shelf Life

Store below 30°C in tightly closed container and away from bright light. Use before expiry date on label.

Reference

- 1) Bauer J.D., Ackermann P.G. and Toro G. (Eds.) 1974, Clinical Laboratory Methods,8th ed. The C.V.Mosby Co.St.Louis.
- 2) Godkar P.B.,1994, Text book of Medical Laboratory Technology, Bhalani Publishing House ,Mumbai, India.

Revision : 1 / 2016



Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory,diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.