



## Vibrio 0129 Differential Disc (10mcg)

DD047

For differentiation of *Vibrio* species based on sensitivity to Vibriostatic agent O129.

### Directions

With a sterile swab, streak the pure, fresh culture of the test organism from sample on a non-selective Blood Agar Plate (containing 0.5% NaCl). Aseptically place both Vibrio 0129 differential disc [10 mcg,(DD047) and 150 mcg (DD048)] on the swabbed plates. Incubate at 35 – 37°C for 24 hours. Observe for zones of inhibition.

### Principle And Interpretation

Shewan and Hodgkiss recognized the sensitivity of *Vibrio* to the vibrio-static agent 0129 (2,4-diamino-6,7-di-isopropyl-pteridine phosphate)(1). 0129 was found to be useful in the differentiation of *Vibrio* from other gram-negative bacteria especially *Aeromonas*, which are characteristically resistant to 0129(2). Even among the genus *Vibrio*, different species show different sensitivities to 0129(3); hence two different concentration discs are to be simultaneously tested to determine the degree of sensitivity of the species. 0129 discs of two concentrations are available: 10-µg and 150-µg. Methods for standardized disc antimicrobial susceptibility testing are employed, with any zone of inhibition around 0129 disks being regarded as sensitive. Medium to be used should be supplemented with 0.5% Sodium Chloride, as sodium ions stimulate the growth of all *Vibrio* species and are required by most.

Interpret the results as follows

Sensitive - Zone of inhibition around both 10µg(DD047)and 150µg(DD048)disc.

Resistant - No zone of inhibition around both the disc (DD047) and (DD048) Partially sensitive - Zone of inhibition around 150µg(DD048)disc and no zone around 10µg(DD047)disc.

### Quality Control

#### Appearance

Vibrio O129 differential filter paper discs of 6 mm diameter containing 10mcg concentration.

#### Cultural response

Cultural response observed on seeded non-selective Blood Agar Plates, with placed Vibrio O129 discs of 10 mcg (DD047) & 150 mcg (DD048),after incubation at 35-37°C for 24 hours.

Organism	Results
<i>Vibrio parahaemolyticus</i> ATCC 17802	Partially sensitive
<i>Aeromonas hydrophila</i> ATCC 7966	Resistant

### Storage and Shelf Life

Store between 2-8°C.For prolonged use store at -20°C. Use before expiry date on label.

### Reference

- 1.Shewan JM, Hodgkiss W. Nature 1954; 63:208-9.
- 2.Isenberg HD, Ed. Clinical microbiology procedures handbook, Vol 1. Washington, DC: ASM, 1994.
- 3.Murray PR, Baron EJ, Pfaller MA, Tenover FC, Tenover FC, Eds. Manual of clinical microbiology. 7th ed. ASM, 1999.
- 4.MacFaddin JF. Biochemical tests for the identification of medical bacteria, 3rd ed. Williams & Wilkins, 2000.

Revision : 1 / 2011

#### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

