



## Neisser's Metachromatic Stains-Kit

K003L

### Intended Use

Neisser's Metachromatic Stains-Kit is used for staining metachromatic granules.

### Composition\*\*

#### Neisser's Methylene Blue (S023)

Ingredients	
Methylene blue	1.000 gm
Ethyl alcohol, 95%	50.000 ml
Glacial acetic acid	50.000 ml
Distilled water	100.000 ml

#### Neutral Red solution (S037)

Ingredients	
Neutral red	1.000 gm
Acetic acid	2.000 ml
Distilled water	1000.000 ml

#### Gram's Iodine (S013)

Ingredients	
Iodine	1.000 gm
Potassium iodide	2.000 gm
Distilled water	300.000 ml

\*\*Formula adjusted, standardized to suit performance parameters

### Directions

- 1) Prepare a smear on clean, dry glass slide.
- 2) Allow it to air dry and fix with gentle heat.
- 3) Stain with Neisser's Methylene Blue (S023) for 3 minutes.
- 4) Wash off with dilute (1:10 in distilled water) Gram's Iodine Solution (S013) and leave some of this solution on the slide for 1 minute.
- 5) Wash with water.
- 6) Counterstain with Neutral Red Solution (S037) for 3 minutes.
- 7) Wash with water and air dry it and observe under oil immersion objective.

### Principle And Interpretation

Well developed granules of volutin (polyphosphate) may be seen in unstained wet preparations as round refractile bodies within the bacterial cytoplasm. With basic dye, they tend to stain more strongly than the rest of the bacterium, and with toluidine blue or methylene blue they stain metachromatically, a reddish purple colour. They are demonstrated most clearly by special methods, such as Albert's and Neisser's, which stain them dark purple but the remainder of the bacterium with a contrasting counterstain. The diphtheria bacillus gives its characteristic volutin-staining reactions best in a young culture (18-24 hours) on a blood or serum medium.

### Type of specimen

Any isolated colony on primary or subculture plates can be isolated from following specimens. Clinical specimen: Blood, urine, CSF, pus, wounds, lesions, body tissues, sputum etc. From environment: Air, water, soil, sludge, waste water, food, dairy samples etc.

### Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines. For food and dairy samples, follow appropriate techniques for sample collection and processing as per guidelines. For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards. Generally the smear is made in laboratory; however, when there is a concern that transport will be delayed or that the preservation for culture will alter the specimen, prepare smear and submit slides to the laboratory.

## Warning and Precautions :

In Vitro diagnostic Use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

## Limitations

1. Overstaining may reduce the contrast between the bacteria and background, or between the cytoplasm and granules.
2. Some strains of *Propionibacterium*, *Actinomyces*, and pleomorphic forms of streptococci may mimic the characteristic stained appearance of *C. diphtheriae*.

## Performance and Evaluation

Performance of the stain is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

## Quality Control

### Microscopic Examination

Metachromatic staining is carried out and staining characteristic of organism is observed under microscope by using oil immersion lens.

### Results

Metachromatic granules : deep blue

The Organisms : pink coloured.

## Storage and Shelf Life

Store between 10 - 30 °C in tightly closed container and away from bright light. Use before expiry date on label. On opening, product should be properly stored in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use.

## Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques.

## Reference

1. Isenberg, H.D. Clinical Microbiology Procedures Handbook. 2nd Edition.
2. Jorgensen., et al. Manual of Clinical Microbiology, American Society for Microbiology, Washington, D.C.
3. Koneman, E.W., et al. Color Atlas and Textbook of Diagnostic Microbiology, J.B. Lippincott Company, Philadelphia, PA.
4. Lillie, R.D. 1977. H.J. Conn's Biological Stains, 9th ed. Williams & Wilkins Company, Baltimore, MD.

Revision : 02 / 2019



In vitro diagnostic medical device



CE Marking



Storage temperature



Do not use if package is damaged



HiMedia Laboratories Pvt. Limited,  
B /4-6 , MIDC, Dindori, Nashik MH  
[www.himedialabs.com](http://www.himedialabs.com)



CE Partner 4U ,Esdoornlaan 13, 3951  
DB Maarn The Netherlands,  
[www.cepartner4u.eu](http://www.cepartner4u.eu)

### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.