



## HI Agar (Heart Infusion Agar)

M169

### Intended Use:

Recommended for isolation and cultivation of a wide variety of fastidious organisms.

### Composition\*\*

Ingredients	Gms / Litre
HM infusion B from #	500.000
Tryptose	10.000
Sodium chloride	5.000
Agar	15.000
Final pH ( at 25°C)	7.4±0.2

\*\*Formula adjusted, standardized to suit performance parameters

# Equivalent to Beef heart, infusion from

### Directions

Suspend 40 grams in 1000 ml of purified / distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. If desired 5% v/v sterile defibrinated blood may be added. Mix well and dispense as desired.

### Principle And Interpretation

Fastidious organisms having exacting nutritional requirement could be cultivated on infusion media, as demonstrated by Huntoon (1). A liquid medium containing an infusion of meat was one of the first media used for the cultivation of bacteria. These infusion media need not be further supplemented by the addition of supplements for cultivation of fastidious bacteria (2). Heart Infusion Agar, containing infusion from HM infusion B is used for the isolation and cultivation of a wide variety of fastidious organisms (3). Heart infusion Agar can also be used for the cultivation of *Vibrio* species (2,4). Heart Infusion Agar can also be supplemented with glucose, horse serum and antibiotics for the cultivation a wide variety of organisms (3). Heart Infusion Agar is used for mass cultivation of organisms for preparation of vaccines. On supplementation of blood, Heart Infusion Agar can be used to study haemolytic reactions (5). This medium was used for isolation and enumeration of haemolytic Streptococci in milk (6).

Tryptose and HM infusion B infusion provide nutritional requirements for the pathogenic bacteria. Sodium chloride maintains the osmotic equilibrium of the medium.

### Type of specimen

Clinical samples - CSF, Blood, Pus, Urine, Feces samples; Dairy samples - Milk samples

### Specimen Collection and Handling

For clinical samples follow appropriate techniques for handling specimens as per established guidelines (2,3). For dairy samples, follow appropriate techniques for sample collection and processing as per guidelines (1,2,8). After use, contaminated materials must be sterilized by autoclaving before discarding.

### Warning and Precautions

In Vitro diagnostic Use only. Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

### Limitations :

1. Further biochemical and serological tests must be carried out for further identification.

## Performance and Evaluation

Performance of the medium is expected as per the direction on the label if it is used as per the direction on the label and stored at recommended temperature!

## Quality Control

### Appearance

Cream to yellow homogeneous free flowing powder

### Gelling

Firm, comparable with 1.5% Agar gel

### Colour and Clarity of prepared medium

Basal medium : Light yellow coloured, clear to slightly opalescent gel After addition of 5-7% w/v sterile defibrinated blood :  
Cherry red coloured, opaque gel forms in Petri plates

### Reaction

Reaction of 4.0% w/v aqueous solution at 25°C. pH : 7.4±0.2

### pH

7.20-7.60

### Cultural Response

Cultural characteristics observed with added 5-7% w/v sterile defibrinated blood, after an incubation at 35-37°C for 18-48 hours.

Organism	Inoculum (CFU)	Growth w/o blood	Recovery w/o blood	Growth with blood	Recovery with blood	Haemolysis
<i>Staphylococcus aureus</i> ATCC 25923	50-100	good-luxuriant	≥70%	luxuriant	≥70%	beta
<i>Neisseria meningitidis</i> ATCC 13090	50-100	luxuriant	≥70%	luxuriant	≥70%	none
<i>Streptococcus pneumoniae</i> ATCC 6303	50-100	good	50-70%	luxuriant	≥70%	alpha
<i>Streptococcus pyogenes</i> ATCC 19615	50-100	good	50-70%	luxuriant	≥70%	beta
<i>Escherichia coli</i> ATCC 25922	50-100	luxuriant	≥70%	luxuriant	≥70%	beta

## Reference

- 1.Huntoon F. M., 1918, J. Inf. Dis., 23:169.
- 2.FDA Bacteriological Analytical Manual, 8th Ed., AOAC International, Gaithersburg, MD.
- 3.Atlas R. M., 2004, Handbook of Microbiological Media, 3rd Ed., CRC Press.
- 4.Downes F. P. and Ito K.,(Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 4th Ed., American Public Health Association, Washington, D.C.
- 5.Murray P. R., Baron J. H., Pfaller M. A., Jorgensen J. H. and Tenover F. C., (Eds.), 2003, Manual of Clinical Microbiology, 8th Ed., American Society for Microbiology, Washington, D.C.
- 6.Diagnostic Procedures and Reagents, 1950, 3rd Edition, 13.

Revision : 1 / 2011

### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.