

# HiGlutaXL™ Minimum Essential Medium Eagle (Alpha modification)

With L-Alanyl-L-Glutamine, Sodium bicarbonate  
Without Deoxyribonucleosides and Ribonucleosides  
1X Liquid Cell Culture Medium

**Product Code: AL081G**

## Product Description :

HiGlutaXL™ medium contains the stabilised dipeptide form of L-glutamine, L-alanyl-L-glutamine. HiGlutaXL™ medium offers several advantages over the conventional glutamine containing media. Dipeptide form prevents the intramolecular cyclization reaction, thus preventing toxic build-up of ammonia. L-alanyl-L-glutamine incorporates L-alanine that protects the alpha amino acid group. Aminopeptidases within the cell break the dipeptide, gradually releasing both L-glutamine and L-alanine for use by the cell. The gradual release of L-glutamine obviates the need to supplement L-glutamine frequently and allows liquid media to be stored at 4°C for longer periods.

AL081G is alpha modification of HiGlutaXL™ Minimum Essential Medium Eagle with L-alanyl-L-glutamine, sodium bicarbonate and non-essential amino acid. It is devoid of deoxyribonucleosides and ribonucleosides. Users are advised to review the literature for recommendations regarding medium supplementation and physiological growth requirements specific for different cell lines.

## Composition :

Ingredients	mg/L
<b>INORGANIC SALTS</b>	
Calcium chloride dihydrate	265.000
Magnesium sulphate anhydrous	97.720
Potassium chloride	400.000
Sodium bicarbonate	2200.000
Sodium chloride	6800.000
Sodium dihydrogen phosphate anhydrous	122.000
<b>AMINO ACIDS</b>	
Glycine	50.000
L-Alanine	25.000
L-Alanyl-L-Glutamine	434.000
L-Arginine hydrochloride	126.000
L-Asparagine monohydrate	44.000
L-Aspartic acid	30.000
L-Cysteine hydrochloride	100.000

L-Cystine dihydrochloride	31.300
L-Glutamic acid	75.000
L-Histidine hydrochloride monohydrate	42.000
L-Isoleucine	52.000
L-Leucine	52.000
L-Lysine hydrochloride	72.500
L-Methionine	15.000
L-Phenylalanine	32.000
L-Proline	40.000
L-Serine	25.000
L-Threonine	48.000
L-Tryptophan	10.000
L-Tyrosine disodium salt	51.900
L-Valine	46.000
<b>VITAMINS</b>	
Choline chloride	1.000
D-Biotin	0.100
D-Ca-Pantothenate	1.000
Folic acid	1.000
L-Ascorbic acid	50.000
Nicotinamide	1.000
Pyridoxal hydrochloride	1.000
Riboflavin	0.100
Thiamine hydrochloride	1.000
Vitamin B12	1.360
i-Inositol	2.000
<b>OTHERS</b>	
D-Glucose	1000.000
Lipoic acid	0.200
Phenol red sodium salt	11.000
Sodium pyruvate	110.000

## Quality Control:

### Appearance

Orangish red colored, clear solution.

### pH

7.00 -7.60

### Osmolality in mOsm/Kg H2O

280.00 -320.00

**Sterility**

No bacterial or fungal growth is observed after 14 days of incubation, as per USP specification.

**Cultural Response**

The growth promotion capacity of the medium is assessed qualitatively by analyzing the cells for the morphology and quantitatively by estimating the cell counts and comparing it with a control medium through minimum three subcultures.

**Endotoxin Content**

NMT 5EU/ml

**Storage and Shelf Life:**

Store at 2-8°C away from bright light.

Shelf life is 18 months.

Use before expiry date given on the product label.

Revision : 1 / 2011

**Disclaimer :**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.