

AS003 Hydrochloric acid abt. 35% pure, Hi-LR™

Product Number	Packing
----------------	---------

AS003	500ML
AS003	2.5L
AS003	5L

Product Information

Product Code	: AS003
Product Name	: Hydrochloric acid abt. 35% pure, Hi-LR™
Molecular Formula	: H C 1
Molecular Weight	: 36.46
CAS No.	: 7647-01-0
EC No.	: 231-595-7
HS Code	: 2806 10 00
Shelf Life	: 3 years

Technical Specification

Appearance	: Clear colourless liquid
Density (d 20/4)	: ~1.19 g/mL
Sulphate (SO ₄)	: <= 0.0005%
Sulphite (SO ₃)	: <= 0.001%
Free chlorine(Cl)	: <= 0.0001%
Heavy metals(as Pb)	: <= 0.0005%
Iron (Fe)	: <= 0.0002%
Arsenic(As)	: <= 0.0001%
Assay (NaOH Titration)	: min. 35.00 %

GHS Safety Information

Hazard Statement(s)	: H314-H335
Precautionary Statement(s)	: P261-P280-P305 + P351 + P338-P310
Signal Word	: Danger
Hazard Pictogram(s)	



Risk and Safety Information

R-Phrase(s)	: 34-37
S-Phrase(s)	: 26-36/37/39-45
WGK	: 1
RTECS	: MW4025000
Storage Temperature(°C)	: Store below 30°C
Hazard Symbol(s)	



Transport Information



Product Information

UN No.	:	1789
Class	:	8
Packaging Group	:	2
Marine Pollutant	:	No
ADR/RID	:	1789 8/PG 2
IMDG	:	1789 8/PG 2
IATA	:	1789 8/PG 2

A-406, Bhaveshwar Plaza, LBS Marg, Mumbai-400 086, India. Phone:+91-022-4095 1919 Fax:+91-022-4095 1920 Gram: STERILITY.
email:info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein