

AS015 Sulfuric acid pure, Hi-LR™

Product Number Packing

AS015	: 500 ML
AS015	: 2.5 L
AS015	: 5 L


Product Information

Product Code	: AS015
Product Name	: Sulfuric acid pure, Hi-LR™
Molecular Formula	: H ₂ SO ₄
Molecular Weight	: 98.08
CAS No.	: 7664-93-9
EC No.	: 231-639-5
HS Code	: 2807 00 00
Shelf Life	: 3 years


Technical Specification

Appearance	: Clear, Colourless Liquid
Density (d 25/4)	: 1.835 - 1.845 g/mL
Chloride (Cl)	: <=0.0005%
Nitrate (NO ₃)	: <=0.0005%
Iron (Fe)	: <=0.0005%
Lead (Pb)	: <= 0.0005%
Residue on ignition	: <= 0.005%
Assay (NaOH Titration)	: min. 95.00%

GHS Safety Information

Hazard Statement(s)	: H314
Precautionary Statement(s)	: P280-P305+P351+P338-P310
Signal Word	: Danger
Hazard Pictogram(s)	

Risk and Safety Information

R-Phrase(s)	: 35
S-Phrase(s)	: 26-30-45
WGK	: 1
Storage Temperature(°C)	: Store below 30°C
Hazard Symbol(s)	

Transport Information

UN No.	: 1830
Class	: 8
Packaging Group	: 2



Product Information

Marine Pollutant	:	No
ADR/RID	:	1830 8/PG 2
IMDG	:	1830 8/PG 2
IATA	:	1830 8/PG 2

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein