

## AS020 Karl Fischer's Reagent, Pyridine free

### Product Number Packing

AS020	500ML
AS020	250ML

### Product Information

Product Code	: AS020
Product Name	: Karl Fischer's Reagent, Pyridine free
Synonym	: Iodine solution methanolic
HS Code	: 3822 00 00
Shelf Life	: 1 years

### Technical Specification

Appearance	: Dark brown liquid
Solubility	: Soluble in water
Determination of the effective factor:	: 1mL = Approx. 5 mg water
Density (d 25/4)	: ~1.15 g/mL

### GHS Safety Information

Hazard Statement(s)	: H225-H301+H311-H331+H360FD+H370
Precautionary Statement(s)	: P201-P210-P233-P240-P280-P302+P352-P304+P340-P403+P235
Signal Word	: Danger
Hazard Pictogram(s)	



### Risk and Safety Information

R-Phrase(s)	: 60-61-11-23/24/25-39/23/24/25
S-Phrase(s)	: 53-16-36/37-45
WGK	: 1
Flash Point(°C)	: 8 °C
Storage Temperature(°C)	: Store below 30°C
Hazard Symbol(s)	



### Transport Information

UN No.	: 1992
Class	: 3
Packaging Group	: II
Marine Pollutant	: No
ADR/RID	: 1992 3/PG II
IMDG	: 1992 3/PG II
IATA	: 1992 3/PG II