

## AS039 Chloroform, Hi-LR™

### Product Number Packing

AS039	: 500ML
AS039	: 2.5L

### Product Information

Product Code	: AS039
Product Name	: Chloroform, Hi-LR™
Synonym	: Trichloromethane, Methylidyne trichloride
Molecular Formula	: CHCl <sub>3</sub>
Molecular Weight	: 119.38
CAS No.	: 67-66-3
EC No.	: 200-663-8
HS Code	: 2903 13 00
Shelf Life	: 3 years

### Technical Specification

Appearance	: Clear colourless liquid
Solubility	: 1 mL miscible with 1 mL of alcohol
FTIR (Liquid film)	: Matches with the standard pattern.
Density (d 25/4)	: 1.487 - 1.489 g/mL
Water (K.F.)	: <= 0.1%
Non-volatile substances	: <= 0.002%
Assay (GC)	: 99.00 - 99.40%

### GHS Safety Information

Hazard Statement(s)	: H302+H332-H315-H319-H336-H351-H361d-H373
Precautionary Statement(s)	: P261-P281-P305+P351+P338
Signal Word	: Warning
Hazard Pictogram(s)	



### Risk and Safety Information

R-Phrase(s)	: 22-38-40-48/20/22
S-Phrase(s)	: 36/37
WGK	: 3
Storage Temperature(°C)	: Store below 30°C
Hazard Symbol(s)	



### Transport Information

UN No.	: 1888
Class	: 6.1
Packaging Group	: 3



## Product Information

Marine Pollutant	:	No
ADR/RID	:	1888 6.1/PG 3
IMDG	:	1888 6.1/PG 3
IATA	:	1888 6.1/PG 3

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919  
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

---

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein