

AS042 Cyclohexane, Hi-LR™

Product Number	Packing
----------------	---------

AS042	500ML
AS042	2.5L

Product Information

Product Code	: AS042
Product Name	: Cyclohexane, Hi-LR™
Molecular Formula	: C ₆ H ₁₂
Molecular Weight	: 84.16
CAS No.	: 110-82-7
EC No.	: 203-806-2
HS Code	: 2902 11 00
EC Index No.	: 601-017-00-1
Shelf Life	: 3 years

Technical Specification

Appearance	: Clear colourless liquid
Density (d 25/4)	: ~0.779 g/mL
Non-volatile substances	: <= 0.005%
Assay (GC)	: min. 99.00 %

GHS Safety Information

Hazard Statement(s)	: H225-H304-H315-H336-H410
Precautionary Statement(s)	: P210-P261-P273-P301+P310-P331-P501
Signal Word	: Danger
Hazard Pictogram(s)	



Risk and Safety Information

R-Phrase(s)	: 11-38-50/53-65-67
S-Phrase(s)	: 9-16-25-33-60-61-62
WGK	: 2
RTECS	: GU6300000
Flash Point(°F)	: -.4 °F
Flash Point(°C)	: -18 °C
Storage Temperature(°C)	: Store below 30°C
Hazard Symbol(s)	



Transport Information

UN No.	: 1145
Class	: 3
Packaging Group	: 2



Product Information

Marine Pollutant	:	YES
ADR/RID	:	1145 3/PG 2
IMDG	:	1145 3/PG 2
IATA	:	1145 3/PG 2

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein