

AS044 Dichloromethane, Hi-LR™

Product Number	Packing
AS044	500ML
AS044	2.5L

Product Information

Product Code	: AS044
Product Name	: Dichloromethane, Hi-LR™
Synonym	: Methylene chloride
Molecular Formula	: $\text{C H}_2 \text{C l}_2$
Molecular Weight	: 84.93
CAS No.	: 75-09-2
EC No.	: 200-838-9
HS Code	: 2903 12 00
EC Index No.	: 602-004-00-3
Shelf Life	: 3 years

Technical Specification

Appearance	: Clear colorless liquid
Density (d 25/4)	: ~1.325 g/mL
Free acid (as HCl)	: $\leq 0.001\%$
Water (K.F.)	: $\leq 0.50\%$
Non-volatile substances	: $\leq 0.005\%$
Assay (GC)	: min.99.0

GHS Safety Information

Hazard Statement(s)	: H351
Precautionary Statement(s)	: P281
Signal Word	: Warning
Hazard Pictogram(s)	



Risk and Safety Information

R-Phrase(s)	: 40
S-Phrase(s)	: 23-24/25-36/37
WGK	: 2
RTECS	: PA8050000
Storage Temperature(°C)	: Store below 30°C
Hazard Symbol(s)	



Transport Information

UN No.	: 1593
Class	: 6.1



Product Information

Packaging Group	:	3
Marine Pollutant	:	YES
ADR/RID	:	1593 6.1/PG 3
IMDG	:	1593 6.1/PG 3
IATA	:	1593 6.1/PG 3

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein