

AS056 Isobutyl methyl ketone, Hi-AR™/ACS

Product Number	Packing
AS056	500ML
AS056	2.5L

Product Information

Product Code	: AS056
Product Name	: Isobutyl methyl ketone, Hi-AR™/ACS
Synonym	: 4-Methyl-2-pentanone; Methyl isobutyl ketone; Isopropylacetone
Molecular Formula	: C ₆ H ₁₂ O
Molecular Weight	: 100.16
CAS No.	: 108-10-1
EC No.	: 203-550-1
HS Code	: 2914 13 00
EC Index No.	: 606-004-00-4
Shelf Life	: 3 years

Technical Specification

Appearance	: Clear colorless liquid
Density (d 25/4)	: ~0.801 g/mL
Refractive index (n 20/D)	: 1.395 - 1.397
Free acid	: <= 0.20 ml N%
Non-volatile substances	: <= 0.005%
Assay (GC)	: min. 98.00%

GHS Safety Information

Hazard Statement(s)	: H225-H319-H332-H335
Precautionary Statement(s)	: P210-P261-P305+P351+P338
Signal Word	: Danger
Hazard Pictogram(s)	



Risk and Safety Information

R-Phrase(s)	: 11-20-36/37-66
S-Phrase(s)	: 9-16-29
WGK	: 1
RTECS	: SA9275000
Flash Point(°F)	: 57.2°F
Flash Point(°C)	: 14°C
Storage Temperature(°C)	: Store below 30°C
Hazard Symbol(s)	



Transport Information

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
 Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such information or with respect to any product, method or apparatus referred to herein



Product Information

UN No.	:	1245
Class	:	3
Packaging Group	:	2
Marine Pollutant	:	No
ADR/RID	:	1245 3/PG 2
IMDG	:	1245 3/PG 2
IATA	:	1245 3/PG 2

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein