

AS064 Petroleum ether, 60 - 80°C, Hi-LR™

Product Number Packing

AS064	: 500ML
AS064	: 2.5L

Product Information

Product Code	: AS064
Product Name	: Petroleum ether, 60 - 80°C, Hi-LR™
CAS No.	: 8032-32-4
EC No.	: 232-453-7
HS Code	: 2710 12 25
Shelf Life	: 3 years

Technical Specification

Appearance	: Clear colorless liquid
Density (d 20/4)	: ~0.68 g/mL
Water (K.F.)	: <= 0.05%
Non-volatile substances	: <= 0.002%

GHS Safety Information

Hazard Statement(s)	: H225-H304-H315-H336-H361-H373-H411
Precautionary Statement(s)	: P210-P261-P273-P281-P301+P310-P331
Signal Word	: Danger
Hazard Pictogram(s)	



Risk and Safety Information

R-Phrase(s)	: 11-38-48/20-51/53-62-65-67
S-Phrase(s)	: 36/37-61-62
WGK	: 2
Flash Point(°F)	: -40 °F
Flash Point(°C)	: -40 °C
Storage Temperature(°C)	: Store below 30°C
Hazard Symbol(s)	



Transport Information

UN No.	: 1268
Class	: 3
Packaging Group	: 1
Marine Pollutant	: Yes
ADR/RID	: 1268 3/PG 1
IMDG	: 1268 3/PG 1
IATA	: 1268 3/PG 1