

AS067 2- Propanol, Hi-LR™**Product Number Packing**

AS067	: 500ML
AS067	: 2.5L
AS067	: 25L


Product Information

Product Code	: AS067
Product Name	: 2- Propanol, Hi-LR™
Synonym	: Isopropyl alcohol; Isopropanol; Dimethyl carbinol
Molecular Formula	: C ₃ H ₈ O
Molecular Weight	: 60.10
CAS No.	: 67-63-0
EC No.	: 200-661-7
HS Code	: 2905 12 00
EC Index No.	: 603-117-00-0
Shelf Life	: 3 years


Technical Specification

Appearance	: Clear colourless liquid
Solubility	: 1 mL miscible in 1 mL of water
FTIR	: Matches with the standard pattern
Refractive index (n 20/D)	: 1.3720 - 1.3820
Density (d 20/4)	: 0.780 - 0.790 g/mL
Non-volatile matter Water (K.F.)	: <= 0.002%
Assay (GC)	: min. 99.00%

GHS Safety Information

Hazard Statement(s)	: H225-H319-H336
Precautionary Statement(s)	: P210-P261-P305+P351+P338
Signal Word	: Danger
Hazard Pictogram(s)	

Risk and Safety Information

R-Phrase(s)	: 11-36-67
S-Phrase(s)	: 7-16-24/25-26
WGK	: 1
RTECS	: NT8050000
Flash Point(°F)	: 53.6 °F
Flash Point(°C)	: 12 °C
Storage Temperature(°C)	: Store below 30°C
Hazard Symbol(s)	



Product Information

Transport Information

Marine Pollutant	:	No
ADR/RID	:	Not Dangerous Goods
IMDG	:	Not Dangerous Goods
IATA	:	Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein