

AS069 Pyridine, Hi-LR™

Product Number Packing

AS069	: 500ML
AS069	: 2.5L

Product Information

Product Code	: AS069
Product Name	: Pyridine, Hi-LR™
Molecular Formula	: C ₅ H ₅ N
Molecular Weight	: 79.10
CAS No.	: 110-86-1
EC No.	: 203-809-9
HS Code	: 2933 31 00
EC Index No.	: 613-002-00-7
Shelf Life	: 3 years

Technical Specification

Appearance	: Clear colourless liquid
Density (d 20/4)	: ~0.98 g/mL
Water (K.F.)	: <= 0.1%
Assay (GC)	: min. 99.00%

GHS Safety Information

Hazard Statement(s)	: H225-H302+H312+H332
Precautionary Statement(s)	: P210-P280
Signal Word	: Danger
Hazard Pictogram(s)	



Risk and Safety Information

R-Phrase(s)	: 11-20/21/22
S-Phrase(s)	: 26-28
WGK	: 2
Flash Point(°F)	: 62.6 °F
Flash Point(°C)	: 17 °C
Storage Temperature(°C)	: Store below 30°C
Hazard Symbol(s)	



Transport Information

UN No.	: 1282
Class	: 3
Packaging Group	: 2
Marine Pollutant	: No



Product Information

ADR/RID : 1282 3/PG 2
IMDG : 1282 3/PG 2
IATA : 1282 3/PG 2

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein