

AS071 Tetrahydrofuran, Hi-AR™/ACS**Product Number Packing**

AS071 : 500ML
AS071 : 2.5L

Product Information

Product Code : AS071
Product Name : Tetrahydrofuran, Hi-AR™/ACS
Molecular Formula : C₄H₈O
Molecular Weight : 72.11
CAS No. : 109-99-9
EC No. : 203-726-8
HS Code : 2932 11 00
EC Index No. : 603-025-00-0
Shelf Life : 3 years

Technical Specification

Appearance : Clear colourless liquid
Density (d 20/4) : ~0.89 g/mL
Peroxides (as H₂O₂) : <= 0.05%
Water (K.F.) : <= 0.05%
Assay (GC) : min. 99.00 %

GHS Safety Information

Hazard Statement(s) : H225-H319-H335-H351
Precautionary Statement(s) : P210-P261-P281-P305+P351+P338
Signal Word : Danger
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 11-19-36/37-40
S-Phrase(s) : 13-16-29-33-36-37-46
WGK : 1
RTECS : LU5950000
Flash Point(°F) : 1.4 °F
Flash Point(°C) : -17 °C
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**

UN No. : 2056
Class : 3



Product Information

Packaging Group	:	2
Marine Pollutant	:	No
ADR/RID	:	2056 3/PG 2
IMDG	:	2056 3/PG 2
IATA	:	2056 3/PG 2

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein