

AS072 Toluene rectified, L.R.**Product Number Packing**

AS072 500ML
AS072 2.5L

Product Information

Product Code : AS072
Product Name : Toluene rectified, L.R.
Synonym : Methylbenzene
Molecular Formula : C₇ H₈
Molecular Weight : 92.14
CAS No. : 108-88-3
EC No. : 203-625-9
HS Code : 2902 30 00
EC Index No. : 601-021-00-3
Shelf Life : 3 years

Technical Specification

Appearance : Clear colourless liquid
Density (d 25/4) : ~0.865 g/mL
Water (K.F.) : <= 0.10%
Non-volatile substances : <= 0.005%
Assay (GC) : min.99.00 %

GHS Safety Information

Hazard Statement(s) : H225-H304-H315-H336-H361d-H373
Precautionary Statement(s) : P210-P261-P281-P301+P310-P331
Signal Word : Danger
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 11-38-48/20-63-65-67
S-Phrase(s) : 36/37-46-62
WGK : 2
RTECS : XS5250000
Flash Point(°F) : 39.2°F
Flash Point(°C) : 4°C
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**

UN No. : 1294



Product Information

Class	:	3
Packaging Group	:	2
Marine Pollutant	:	No
ADR/RID	:	1294 3/PG 2
IMDG	:	1294 3/PG 2
IATA	:	1294 3/PG 2

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein