

AS073 Toluene low in sulfur, L.R.

Product Number	Packing
----------------	---------

AS073	500ML
AS073	2.5L

Product Information

Product Code	: AS073
Product Name	: Toluene low in sulfur, L.R
Synonym	: Methylbenzene
Molecular Formula	: C ₇ H ₈
Molecular Weight	: 92.14
CAS No.	: 108-88-3
EC No.	: 203-625-9
HS Code	: 2902 30 00
EC Index No.	: 601-021-00-3
Shelf Life	: 3 years

Technical Specification

Appearance	: Clear colourless liquid
Density (d 25/4)	: ~0.865 g/mL
Sulfur compounds (C2S)	: <= 0.0005%
Water (K.F.)	: <= 0.04%
Non-volatile substances	: <= 0.003%
Assay (GC)	: min. 99.00 %

GHS Safety Information

Hazard Statement(s)	: H225-H304-H315-H336-H361d-H373
Precautionary Statement(s)	: P210-P261-P281-P301+P310-P331
Signal Word	: Danger
Hazard Pictogram(s)	



Risk and Safety Information

R-Phrase(s)	: 11-38-48/20-63-65-67
S-Phrase(s)	: 36/37-46-62
WGK	: 2
RTECS	: XS5250000
Flash Point(°F)	: 39.2°F
Flash Point(°C)	: 4°C
Storage Temperature(°C)	: Store below 30°C
Hazard Symbol(s)	



Transport Information



Product Information

UN No.	:	1294
Class	:	3
Packaging Group	:	2
Marine Pollutant	:	No
ADR/RID	:	1294 3/PG 2
IMDG	:	1294 3/PG 2
IATA	:	1294 3/PG 2

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein