

## AS075 Trichloroethylene, L.R.

### Product Number Packing

AS075	500ML
AS075	2.5L

### Product Information

Product Code	: AS075-500ML
Product Name	: Trichloroethylene, L.R.
Synonym	: Trichloroethene; Ethylene trichloride; 1,1,2-Trichloroethene
Molecular Formula	: $C_2HCl_3$
Molecular Weight	: 131.39
CAS No.	: 79-01-6
EC No.	: 201-167-4
HS Code	: 2903 22 00
EC Index No.	: 602-027-00-9
Shelf Life	: 3 years


### Technical Specification

Appearance	: Clear colourless liquid
Density (d 25/4)	: ~1.463 g/mL
Water (K.F.)	: <= 0.05%
Non-volatile substances	: <= 0.004%
Assay (GC)	: min. 99.00 %

### GHS Safety Information

Hazard Statement(s)	: H315-H319-H336-H341-H350-H412
Precautionary Statement(s)	: P201-P261-P273-P281-P305+P351+P338-P308+P313
Signal Word	: Danger
Hazard Pictogram(s)	

### Risk and Safety Information

R-Phrase(s)	: 45-36/38-52/53-67-68
S-Phrase(s)	: 53-45-61
WGK	: 3
RTECS	: KX4550000
Storage Temperature(°C)	: Store below 30°C
Hazard Symbol(s)	

### Transport Information

UN No.	: 1710
Class	: 6.1
Packaging Group	: 3



## Product Information

Marine Pollutant	:	No
ADR/RID	:	1710 6.1/PG 3
IMDG	:	1710 6.1/PG 3
IATA	:	1710 6.1/PG 3

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919  
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

---

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein