

AS105 Carbinol, Hi-LR™

Product Number	Packing
----------------	---------

AS105	500ML
AS105	2.5L

Product Information

Product Code	: AS105
Product Name	: Carbinol, Hi-LR™
Molecular Formula	: CH ₃ OH
Molecular Weight	: 32.04
CAS No.	: 67-56-1
EC No.	: 200-659-6
HS Code	: 2905 11 00
Shelf Life	: 3 years

Technical Specification

Appearance	: Clear colourless liquid
Solubility	: 1 mL miscible in 1 mL of water
Density (d 25/4)	: ~0.791 g/mL
Non volatile substances	: <= 0.005%
Water (K.F.)	: <= 0.20%
Aldehydes and ketones [(CH ₃) ₂ CO]	: <= 0.10%
Assay (GC)	: min. 99.00 %


GHS Safety Information

Hazard Statement(s)	: H225-H301-H311-H331-H370
Precautionary Statement(s)	: P210-P260-P280-P301+P310-P311
Signal Word	: Danger
Hazard Pictogram(s)	



Risk and Safety Information

R-Phrase(s)	: 11-23/24/25-39/23/24/25
S-Phrase(s)	: 7-16-36/37-45
WGK	: 1
RTECS	: PC1400000
Flash Point(°F)	: 51.8 °F
Flash Point(°C)	: 11 °C
Storage Temperature(°C)	: Store below 30°C
Hazard Symbol(s)	



Transport Information



Product Information

UN No.	:	1230
Class	:	3
Packaging Group	:	2
Marine Pollutant	:	No
ADR/RID	:	1230 3/PG 2
IMDG	:	1230 3/PG 2
IATA	:	1230 3/PG 2

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein