

**AS109 Isoamyl alcohol, Hi-AR™/ACS****Product Number Packing**

AS109 : 2.5L  
AS109 : 500ML

**Product Information**

Product Code : AS109  
Product Name : Isoamyl alcohol, Hi-AR™/ACS  
Molecular Formula : C<sub>5</sub>H<sub>12</sub>O  
Molecular Weight : 88.15  
CAS No. : 123-51-3  
EC No. : 204-633-5  
EC Index No. : 603-006-00-7  
Shelf Life : 3 years

**Technical Specification**

Appearance : Clear, Colourless Liquid.  
Identity (IR) : Passes test  
Density (d 20/4) : 0.806 - 0.809 g/mL  
Non-volatile substances : <= 0.01%  
Assay (GC) : min. 98.00 %

**GHS Safety Information**

Hazard Statement(s) : H226-H332-H335  
Precautionary Statement(s) : P261  
Signal Word : Warning  
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 10-20-37-66  
S-Phrase(s) : 46  
WGK : 1  
Flash Point(°F) : 109.4 °F  
Flash Point(°C) : 43 °C  
Storage Temperature(°C) : Store below 30°C  
Hazard Symbol(s)

**Transport Information**

UN No. : 1105  
Class : 3  
Packaging Group : 3  
Marine Pollutant : No



## Product Information

ADR/RID : 1105 3/PG 3  
IMDG : 1105 3/PG 3  
IATA : 1105 3/PG 3

A-406, Bhaveshwar Plaza, LBS Marg, Mumbai-400 086, India. Phone:+91-022-4095 1919 Fax:+91-022-4095 1920 Gram: STERILITY.  
email:info@himedialabs.com. Visit us at our website: [www.himedialabs.com](http://www.himedialabs.com)

---

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein