

AS121 Dimethyl sulphoxide, Hi-LR™

Product Number Packing

AS121 : 500ML

Product Information

Product Code : AS121
Product Name : Dimethyl sulphoxide, Hi-LR™
Molecular Formula : C₂H₆OS
Molecular Weight : 78.13
CAS No. : 67-68-5
EC No. : 200-664-3
HS Code : 2930 90 99
Shelf Life : 3 years

Technical Specification

Appearance : Clear colorless liquid
Refractive index (n 20/D) : 1.475 - 1.480
Density (d 20/4) : 1.099 - 1.101
Acidity (as CH₃COOH) : <= 0.001%
Water (K.F.) : <= 0.2%
Assay (GC) : min. 99.00%

Risk and Safety Information

WGK : 1
RTECS : PV6210000
Flash Point(°F) : 188.6 °F
Flash Point(°C) : 87 °C
Storage Temperature(°C) : Store below 30°C

Transport Information

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods