



Product Information

CG272 15ml Centrifuge Tube, with RED screw cap

Product Information

Product Name : 15ml Centrifuge Tube, with RED screw cap

Technical Specification

Description : Centrifuge tube having conical bottom with red screw cap, printed graduations and made up of polypropylene. Autoclavable.

Dimensions : 16 mm x 120 mm

Volume capacity : 15 ml

Temperature tolerance test : Withstands steam sterilization at 121°C for 30 minutes and not to be used in dry heat sterilization.

Use : For centrifugation.

Transport Information

ADR/RID : Not Dangerous Goods

IMDG : Not Dangerous Goods

IATA : Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such information or with respect to any product, method or apparatus referred to herein