

CG279

2 mL Moulded Graduation, Centrifuge Tubes with Lid

Product Number

Packing

CG279 : 1 x 1

Product Information

Product Code : CG279
Product Name : 2 mL Moulded Graduation, Centrifuge Tubes with Lid

Technical Specification

Description : Autoclavable round bottom centrifuge tube with lid having moulded graduation made up of polypropylene material.
Clarity : Clear without any scratches and fibers.
Capacity : 2 mL
Product alignment : Lid closer of the tube is tight.
Temperature tolerance test : Withstand steam sterilization at 121°C for 30 minutes and not to be used in dry heat sterilization.
Use : For centrifugation.

Transport Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods