



# Technical Data

## Campylobacter Supplement - II (Butzler)

FD007

An antibiotic supplement recommended for the selective isolation of *Campylobacter* species.

### Composition

Per vial sufficient for 500 ml medium

*Ingredients	Concentration
Bacitracin	12500IU
Colistin sulphate	5000Units
Cycloheximide	25mg
Cephazolin sodium	7.500mg
Novobiocin	2.500mg

### Directions:

Rehydrate the contents of one vial aseptically with 5 ml of 50% aqueous ethanol and mix well to dissolve. Avoid frothing of the solution. Aseptically add the contents to 500 ml of sterile, molten, cooled (45-50°C) Columbia Blood Agar Base [M144](#) / Columbia Blood Agar Base, HiVeg™ [MV144](#) / Columbia Blood Agar Base w/ 1% Agar [M144A](#) / Columbia Blood Agar Base w/ 1% Agar, HiVeg™ [MV144A](#) / Blood Agar Base No. 2 [M834](#) / Blood Agar Base No. 2 w/ 1.2% Agar [M834A](#) / Blood Agar Base No. 2, HiVeg™ [MV834](#) / Blood Agar Base No. 2 w/ 1.2% Agar, HiVeg™ [MV834A](#) / Brucella Agar Base [M074](#) / Brucella HiVeg™ Agar Base [MV074](#) along with 5% v/v defibrinated Sheep blood or 5-7% v/v laked horse blood and rehydrated contents of 1 vial of Campylobacter Growth Supplement [FD009](#) . For 470 ml of Modified Brucella Agar Base [M1268](#) , add 1 vial of rehydrated contents of FD007 along with 5% v/v defibrinated Sheep blood. Mix gently and pour into sterile petri plates.

### Storage and Shelf Life

Store at 2 - 8°C. Use before expiry date on label.

Revision : 1 / 2012

\* Not For Medicinal Use

### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.