

Technical Data

Campylobacter Supplement- III (Skirrow)

FD008

An antibiotic supplement recommended for the selective isolation of *Campylobacter* species.

Composition

Per vial sufficient for 500 ml medium

*Ingredients Concentration
Vancomycin 5mg
Polymyxin B sulphate 1250IU
Trimethoprim 2.500mg

Directions:

Rehydrate the contents of one vial aseptically with 2 ml sterile distilled water and mix well to dissolve. Avoid frothing of the solution. Aseptically add the rehydrated contents to 500 ml of sterile, molten, cooled (45-50°C) Columbia Blood Agar Base M144 / Columbia Blood Agar Base, HiVegTM MV144 / Columbia Blood Agar Base w/ 1% Agar M144A / Columbia Blood Agar Base w/ 1% Agar, HiVegTM MV144A / Blood Agar Base No. 2 M834 / Blood Agar Base No. 2, HiVegTM / Blood Agar Base No. 2 w/ 2% Agar M834A / Blood Agar Base No. 2 w/ 2% Agar, HiVegTM MV834A Campylobacter Agar Base M994 / Campylobacter HiVegTM Agar Base MV994 with 5-7% v/v lysed defibrinated horse blood/ Brucella MV074 / Brucella Agar Base, Modified Agar Base M074 / Brucella HiVegTM Agar Base M074A / Brucella HiVegTM Agar Base, Modified MV074A with 5% inactivated horse serum. Mix gently and pour into sterile petri plates.

Storage and Shelf Life

Store at 2-8°C. Use before the expiry date on the label.

Revision: 1 / 2012

* Not For Medicinal Use

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia[™] publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia[™] Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.