



Technical Data

GC Supplement w/ Antibiotics

FD021

An antibiotic and enrichment supplement recommended for the selective isolation of pathogenic *Neisseria* .

Composition

Per vial sufficient for 500 ml medium

*Ingredients	Concentration
Yeast autolysate	5g
Colistin methane sulphionate	3.750mg
Dextrose	0.750g
Trimethoprim	2.500mg
Sodium bicarbonate	0.075g
Nystatin	6250Units
Vancomycin	1.500mg

Directions:

Rehydrate the contents of 1 vial aseptically with 15 ml of sterile distilled water. Mix well and add aseptically to 500 ml of sterile, molten, cooled (45-50°C) GC Agar Base [M434](#) / GC HiVeg™ Agar Base [MV434](#) / Thayer Martin Medium Base [M413](#) / Thayer Martin HiVeg™ Medium Base [MV413](#) along with separately prepared Hemoglobin [FD022](#). Mix well and pour into sterile petri plates.

Storage and Shelf Life

Store at 2-8°C. Use before expiry date on the label.

Revision : 1 / 2012

* Not For Medicinal Use

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.