



Technical Data

Linco T Supplement

FD026

An antibiotic supplement recommended for the selective isolation and cultivation of *Neisseria* species.

Composition

Per vial sufficient for 500 ml medium

*Ingredients	Concentration
Lincomycin	0.500mg
Colistin	3mg
Amphotericin B	0.500mg
Trimethoprim lactate	3.250mg

Directions:

Rehydrate the contents of 1 vial aseptically with 10 ml of sterile distilled water. Mix gently, avoiding frothing to dissolve and aseptically add it to 240 ml of sterile, molten, cooled (45-50°C), GC Agar Base [M434](#) / GC HiVeg™ Agar Base [MV434](#) along with 250 ml of sterile 2% w/v Haemoglobin solution [FD022](#) . Also add sterile rehydrated contents of 1 vial of Yeast Autolysate Supplement [FD027](#) . Mix gently and pour into sterile petri plates.

Storage and Shelf Life

Store at 2-8°C. Use before expiry date on the label.

Revision : 1 / 2012

* Not For Medicinal Use

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.