



# Technical Data

## Amylose Azure (1.5 gm per vial)

FD060

This supplement is recommended for the cultivation, differentiation and enumeration of Enterococcus species based on starch hydrolysis and production of fluorescence.

### Composition

Per vial sufficient for 1000 ml medium

#### Ingredients

Amylose azure

#### Concentration

1.500g

### Directions:

Rehydrate the contents of two vials aseptically with 10 ml of sterile distilled water. Mix well and aseptically add it to sterile, molten, cooled (45-50°C) 970 ml of FGTC Agar Base [M1171](#) / FGTC HiVeg™ Agar Base [MV1171](#) along with rehydrated contents of 1 vial each of FGTC Antibiotic Supplement [FD079](#) and Sodium Bicarbonate Solution [FD080](#). Mix well before pouring into sterile petri plates.

### Storage and Shelf Life

Store between 10-30°C. Use before expiry date on the label.

Revision : 1 / 2012

#### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.