



Technical Data

Willis and Hobb's Supplement

FD156

An antibiotic supplement for isolation of *Clostridium perfringens* and *Clostridium botulinum* species from food.

Composition

Per vial sufficient for 475 ml medium

*Ingredients

Neomycin sulphate

Concentration

33mg

Directions:

Rehydrate the contents of one vial aseptically with 5.0 ml of sterile distilled water and aseptically add it to 400 ml of sterile, molten, cooled (45-50°C) Willis and Hobbs Medium Base [M1375](#) along with 15 ml of Egg Yolk Emulsion [FD045](#) and 60 ml of sterile skim milk.

Storage and Shelf Life

Store at 2-8°C. Use before expiry date on the label.

Revision : 1 / 2012

* Not For Medicinal Use

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.