

## Liver Meat Agar, Granulated

GM1001

Liver Meat Agar, granulated is used for cultivation of fastidious anaerobic microorganisms.

### Composition\*\*

Ingredients	Gms / Litre
Meat liver infusion powder	20.000
Dextrose	0.750
Starch	0.750
Sodium sulphite	1.200
Ferric ammonium citrate	0.500
Agar	11.000
Final pH ( at 25°C)	7.6±0.2

\*\*Formula adjusted, standardized to suit performance parameters

### Directions

Suspend 34.2 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well and pour into sterile Petri plates.

### Principle And Interpretation

Anaerobic bacteria live in an oxygen-free environment. Some anaerobic bacteria actually die if oxygen is present, while others fail to grow and multiply (1). Meat liver base provides adequate degree of anaerobiosis and is also rich source of growth nutrients, which enables even the strict and fastidious anaerobes to grow well.

Some anaerobes (e.g. certain *Clostridium* species) reduce the sulphite present in the medium to hydrogen sulphide (H<sub>2</sub>S) which is indicated by the blackening of colonies due to presence of ferric ammonium citrate. Inoculation can be performed by the pour plate method or by surface smearing

### Quality Control

#### Appearance

Light yellow to light brown coloured granular medium

#### Gelling

Firm, comparable with 1.1% Agar gel

#### Colour and Clarity of prepared medium

Brown coloured opalescent gel with suspended particles forms in Petri plates.

#### Reaction

Reaction of 3.42% w/v aqueous solution at 25°C. pH : 7.6±0.2

#### pH

7.40-7.80

#### Cultural Response

Cultural characteristics observed under anaerobic condition, after an incubation at 35-37°C for 18-48 hours.

Organism	Inoculum (CFU)	Growth	Recovery	H <sub>2</sub> S
<b>Cultural Response</b>				
<i>Clostridium perfringens</i> ATCC 12924	50-100	luxuriant	≥50%	positive
<i>Clostridium tetani</i> ATCC 10779	50-100	luxuriant	≥50%	positive
<i>Escherichia coli</i> ATCC 25922	50-100	luxuriant	≥50%	negative
<i>Proteus mirabilis</i> ATCC 25933	50-100	luxuriant	≥50%	negative or weakly positive

---

<i>Clostridium botulinum</i> ATCC 25763	50-100	luxuriant	>=50%	positive
<i>Bacteroides vulgatus</i> ATCC 8482	50-100	good-luxuriant	>=50%	negative

### Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label.

### Reference

1. Alcamo E. I., 2001, Fundamentals of Microbiology, 6th Ed., Jones and Bartlett Publishers

Revision : 00 / 2014



#### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.