

GRM1414 Silver metal powder**Product Number Packing**

GRM1414 : 25G

Product Information

Product Code : GRM1414
Product Name : Silver metal powder
Atomic Formula : Ag
Atomic Weight : 107.87
CAS No. : 7440-22-4
EC No. : 231-131-3
HS Code : 7106 10 00
Shelf Life : 4 years

Technical Specification

Appearance : Silver to silver grey lustrous powder or granules
Solubility : 33.3 mg soluble in 1 mL of nitric acid
Assay (AT) : min. 99.99 %

Risk and Safety Information

WGK : 3
RTECS : VW3500000
Storage Temperature(°C) : Store below 30°C

Transport Information

Marine Pollutant : -
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods