

GRM3051 Manganese metal powder

Product Number Packing

GRM3051 : 250G

Product Information

Product Code : GRM3051
Product Name : Manganese metal powder
CAS No. : 7439-96-5
Shelf Life : 4 years

Technical Specification

Appearance : Steel grey, lustrous powder.
Solubility : Insoluble in water.
Assay (EDTA Titration) : min. 99.00 %

Risk and Safety Information

R-Phrase(s) : 11-15
S-Phrase(s) : 43
RTECS : OO9275000
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)



Transport Information

UN No. : 3208
Class : 4.3
Packaging Group : 1
ADR/RID : 3208 4.3/PG 1
IMDG : 3208 4.3/PG 1
IATA : 3208 4.3/PG 1