

GRM3313 Cobalt metal powder**Product Number Packing**

GRM3313	: 100G
GRM3313	: 250G


Product Information

Product Code	: GRM3313
Product Name	: Cobalt metal powder
Atomic Formula	: Co
Atomic Weight	: 58.93
CAS No.	: 7440-48-4
EC No.	: 231-158-0
HS Code	: 8105 20 00
EC Index No.	: 027-001-00-9
Shelf Life	: 4 years


Technical Specification

Appearance	: Orange to dark grey to black powder
Solubility	: 33.3 mg soluble in 1 mL of dilute nitric acid
Assay (EDTA Titration)	: min. 99.50%

GHS Safety Information

Hazard Statement(s)	: H317-H334-H413
Precautionary Statement(s)	: P261-P280-P342+P311
Signal Word	: Danger
Hazard Pictogram(s)	

Risk and Safety Information

R-Phrase(s)	: 42/43-53
S-Phrase(s)	: 22-24-37-61
WGK	: nwg
RTECS	: GF8750000
Storage Temperature(°C)	: Store below 30°C
Hazard Symbol(s)	

Transport Information

UN No.	: 3089
Class	: 4.1
Packaging Group	: 2
Marine Pollutant	: No



Product Information

ADR/RID : 3089 4.1/PG2
IMDG : 3089 4.1/PG2
IATA : 3089 4.1/PG2

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein