

GRM4504 Tellurium metal powder**Product Number Packing**

GRM4504 : 100G

Product Information

Product Code : GRM4504
Product Name : Tellurium metal powder
Atomic Formula : Te
Atomic Weight : 127.60
CAS No. : 13494-80-9
EC No. : 236-813-4
HS Code : 2804 50 90
Shelf Life : 4 years

Technical Specification

Appearance : Grey to black powder
Solubility : 33.3 mg soluble in 1 mL of sulphuric acid (Opaque solution)
ICP : Confirms Tellurium metal
Assay (Based on Trace metal analysis): min. 99.90%

GHS Safety Information

Hazard Statement(s) : H301
Precautionary Statement(s) : P301+P310
Signal Word : Danger
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 25
S-Phrase(s) : 45
WGK : 3
RTECS : WY2625000
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**

Marine Pollutant : No
ADR/RID : Not dangerous goods
IMDG : Not dangerous goods
IATA : Not dangerous goods