

GRM7182 Lithium metal (in liquid paraffin)**Product Number Packing**

GRM7182 : 500G

Product Information

Product Code : GRM7182
Product Name : Lithium metal (in liquid paraffin)
Atomic Formula : Li
Atomic Weight : 6.94
CAS No. : 7439-93-2
EC No. : 231-102-5
HS Code : 2805 19 90
EC Index No. : 003-001-00-4
Shelf Life : 4 years

Technical Specification


Appearance : Silvery white to grey metal (in liquid paraffin)
Solubility : Reacts violently with water
Assay (ICP) : min. 99.00%

GHS Safety Information

Hazard Statement(s) : H260-H314
Precautionary Statement(s) : P223-P231+P232-P280-P305+P351+P338-P370+P378-P422
Signal Word : Danger
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 14/15-34
S-Phrase(s) : 8-43-45
WGK : 2
RTECS : OJ5540000
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**

UN No. : 1415
Class : 4.3
Packaging Group : 1
Marine Pollutant : No
ADR/RID : 1415 4.3/PG 1
IMDG : 1415 4.3/PG 1



Product Information

IATA

: 1415 4.3/PG 1

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein