

GRM740 Aluminium metal fine powder

Product Number Packing

GRM740 : 500G

Product Information

Product Code : GRM740
Product Name : Aluminium metal fine powder
Atomic Formula : Al
Atomic Weight : 26.98
CAS No. : 7429-90-5
EC No. : 231-072-3
HS Code : 7603 10 00
Shelf Life : 4 years

Technical Specification

Appearance : Tin-white to grey powder
Solubility : Reacts with dilute hydrochloric acid with the evolution of hydrogen
Assay (EDTA Titration) : 99.00 - 101.00 %

GHS Safety Information

Hazard Statement(s) : H250-H261
Precautionary Statement(s) : P222-P231+P232-P422
Signal Word : Danger
Hazard Pictogram(s)



Risk and Safety Information

R-Phrase(s) : 15 -17
S-Phrase(s) : 43-7/8
WGK : nwg
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)



Transport Information

UN No. : 1396
Class : 4.3
Packaging Group : 2
EMS Code : F-G, S-O
Marine Pollutant : No
ADR/RID : 1396 4.3/PG 2



Product Information

IMDG : 1396 4.3/PG 2
IATA : 1396 4.3/PG 2

A-406, Bhaveshwar Plaza, LBS Marg, Mumbai-400 086, India. Phone:+91-022-4095 1919 Fax:+91-022-4095 1920 Gram: STERILITY.
email:info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein