

**GRM7563 Tin metal powder, Pure****Product Number Packing**

GRM7563 : 100G

**Product Information**

Product Code : GRM7563  
Product Name : Tin metal powder, Pure  
Atomic formula : Sn  
Atomic Weight : 118.71  
CAS No. : 7440-31-5  
EC No. : 231-141-8  
HS Code : 8007 00 80  
Shelf Life : 4 years

**Technical Specification**

Appearance : Grey powder  
Solubility : Insoluble in water  
Particle size : min. 80.00% passing through 325 mesh (44 µm)  
Assay (EDTA Titration) : 99.00 - 102.00%

**GHS Safety Information**

Hazard Statement(s) : H319-H335  
Precautionary Statement(s) : P261-P305+P351+P338  
Signal Word : Warning  
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 36/37  
S-Phrase(s) : 26  
WGK : 3  
RTECS : XP7320000  
Storage Temperature(°C) : Store below 30°C  
Hazard Symbol(s)

**Transport Information**

Marine Pollutant : No  
ADR/RID : Not Dangerous Goods  
IMDG : Not Dangerous Goods  
IATA : Not Dangerous Goods