

GRM7769 Aluminium metal pellets

Product Number Packing

GRM7769 : 500G

Product Information

Product Code : GRM7769
Product Name : Aluminium metal pellets
Atomic Formula : Al
Atomic Weight : 26.98
CAS No. : 7429-90-5
EC No. : 231-072-3
HS Code : 7601 10 00

Technical Specification

Appearance : Tin white to grey metal pellets
Solubility : 33.3 mg soluble in 1 mL of dilute hydrochloric acid
Assay (EDTA Titration) : min. 99.00 %

GHS Safety Information

Hazard Statement(s) : H400
Precautionary Statement(s) : P273
Signal Word : Warning
Hazard Pictogram(s)



Risk and Safety Information

WGK : 3
RTECS : BD0330000
Storage Temperature(°C) : Store below 30°C

Transport Information

Marine Pollutant : Yes
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods