

MB094 Agarose, For Pulse Field Gel Electrophoresis (PFGE), For Molecular Biology**Product Number****Packing**

MB094	: 100G
MB094	: 1KG
MB094	: 500G

Product Information

Product Code	: MB094
Product Name	: Agarose, For Pulse Field Gel Electrophoresis (PFGE), For Molecular Biology
CAS No.	: 9012-36-6
EC No.	: 232-731-8
HS Code	: 3913 90 00
Other Information	: 1) Gel point : The temperature at which an aqueous agarose solution forms a gel as it cools. 2) Formulated for fast resolution of large DNA by pulsed-field gel electrophoresis (PFGE) due to it's higher exclusion limits.
Shelf Life	: 4 years

Technical Specification

Appearance	: White to off-white powder
Solubility	: 33.3 mg soluble in 1 mL of hot water
DNases	: None detected
RNases	: None detected
Gel temperature (10% Gel)	: 41 - 43°C.
Water (K.F.)	: <=10.0%
EEO	: 0.04 - 0.10

Risk and Safety Information

WGK	: 3
Storage Temperature(°C)	: Store below 30°C

Transport Information

Marine Pollutant	: No
ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods