

MB175 Agarose special, Medium EEO, For Molecular Biology

Product Number Packing

MB175	: 10G
MB175	: 100G
MB175	: 1KG

Product Information

Product Code	: MB175
Product Name	: Agarose special, Medium EEO, For Molecular Biology
CAS No.	: 9012-36-6
EC No.	: 232-731-8
HS Code	: 3913 90 00

Technical Specification

Appearance	: White powder
Solubility	: 33.3 mg soluble in 1 mL of hot water
Residue on ignition (ash)	: <=0.5%
Sulphate (SO ₄)	: <=0.5%
Gel point (1.5%)	: 34.5 - 37.5°C
EEO	: 0.15 - 0.2
DNases	: None detected
RNases	: None detected

Risk and Safety Information

WGK	: 3
Storage Temperature(°C)	: Store below 30°C

Transport Information

Marine Pollutant	: No
ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods