

Centrifuge Tube, 50 ml

PW1309

These conical bottom centrifuge tubes are autoclavable having orange cap. Material is polypropylene. Handy packaging makes our centrifuge tubes the economical choice.

Application : Used in the laboratories like Molecular biology laboratory, microbiology laboratory, Tissue culture work, Biochemistry and cell biology in centrifugation procedures.

| Product Name | Product Code | Description | Capacity (ml) | Maximum RCF value |
|-----------------|--------------|--|---------------|-------------------|
| Centrifuge Tube | PW1309 | Sterile, Autoclavable , Conical Bottom with Orange cap and printed graduations, large writing area. DNase / RNase free. Pyrogen free | 50 .00 | 20000 g |

Product Features:

- Conical Bottom with Orange Cap.
- Sterile, Autoclavable.
- Printed Graduation

Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.