



## Malt Extract Powder, Refined

RM004B

It is rich in carbohydrates and is ideally suitable for use in media employed for cultivation of fungi such as, Malt Agar, Malt Extract Agar, Yeast Malt Agar, etc.

### Principle And Interpretation

Malt Extract Powder is prepared from aqueous extract of sprouted malt grains and dried at low temperature to preserve nutrients present in the form of carbohydrates and nitrogenous substances.

### Quality Control

#### Appearance

Cream to brownish yellow Homogenous Free flowing powder

#### Solubility

Freely soluble in distilled/purified water, insoluble in alcohol.

#### Clarity

1% w/v aqueous solution remains clear without haziness after autoclaving at 15 lbs pressure (121°C) for 15 minutes.

#### Reaction

Reaction of 10% w/v aqueous solution at 25°C.

#### pH

5.00- 7.00

#### Microbial Load:

##### Total aerobic microbial count (cfu/gm)

By plate method when incubated at 30-35°C for not less than 3 days.

Bacterial Count : <= 2000 CFU/gram

##### Total Yeast and mould count (cfu/gm)

By plate method when incubated at 20-25°C for not less than 5 days.

Yeast & mould Count : <= 100 CFU/gram

#### Test for Pathogens

1. E.coli-Negative in 10 gms of sample  
2. Salmonella species-Negative in 10 gms of sample  
3. Pseudomonas aeruginosa-Negative in 10 gms of sample  
4. Staphylococcus aureus- Negative in 10 gms of sample  
5. C.albicans- Negative in 10 gms of sample  
6. Clostridia- Negative in 10 gms of sample

#### Cultural response

Cultural response observed after an incubation at 25-30°C for 48-72 hours by preparing Malt Extract Agar (M137) using Malt Extract as an ingredient.

#### Cultural Response

Organism	Growth
<b>Cultural response</b>	
<i>Candida albicans</i> ATCC 10231	Luxuriant
<i>Saccharomyces cerevisiae</i> ATCC 9763	Luxuriant
<i>Aspergillus brasiliensis</i> ATCC 16404	Luxuriant

**Chemical Analysis**

Protein Content	$\geq 5.0\%$
Carbohydrate Content	$\geq 70\%$
Sodium chloride	$\leq 1.0\%$
Loss on drying	$\leq 6.0\%$
Residue on ignition	$\leq 4.50\%$

**Storage and Shelf Life**

Store below 30°C. Use before expiry date on the label.

**Disclaimer :**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.