

RM1352 Calcium metal**Product Number Packing**

RM1352 : 100G
RM1352 : 5X100G

Product Information

Product Code : RM1352
Product Name : Calcium metal
Atomic Formula : Ca
Atomic Weight : 40.08
CAS No. : 7440-70-2
EC No. : 231-179-5
HS Code : 2805 12 00
Shelf Life : 4 years

Technical Specification

Appearance : Silver or grey or gold granules
Solubility : 33.3 mg soluble in 1 mL of conc. hydrochloric acid
Assay (EDTA Titration) : 98.50 - 101.50%

GHS Safety Information

Hazard Statement(s) : H261
Precautionary Statement(s) : P231+P232-P422
Signal Word : Danger
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 15
S-Phrase(s) : 8-24/25-43
WGK : 1
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**

UN No. : 1401
Class : 4.3
Packaging Group : 2
Marine Pollutant : No
ADR/RID : 1401 4.3/PG 2
IMDG : 1401 4.3/PG 2
IATA : 1401 4.3/PG 2