
Copper metal foil, Hi-LR™, ~ 0.1 mm thick

RM3044

Product Identifier

CAS No.	:	7440-50-8
EC No.	:	231-159-6
Atomic Formula	:	Cu
Atomic Weight	:	63.55
HS Code	:	7410 11 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Reddish lustrous metal foil of about 0.1 mm thickness
Solubility	:	33.3 mg soluble in 1 mL of nitric acid
Nitric acid-insoluble matter	:	<= 0.01%
Assay (Iodometry/EDTA Titration)	:	99.00 - 102.00%

Safety Information

UN No.	:	3077
Class	:	9
Packing Group	:	III
RTECS	:	GL5325000
WGK	:	3